

# Product Catalogue



# incidental



MADE IN  
TÜRKİYE

**ENDOART®**  
FINAL TOUCH

**ruby | dent**





**ACTION**

SX	S1	S2	F1	F2	F3	F4	F5		SX	S1	S2	F1	F2	F3	F4	F5	
Starter		Shaping		Finishing					Starter		Shaping		Finishing				
Torque: 3.0 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm

↑ Excellent Cutting Efficiency     ↻ Maximum File Flexibility

EndoArt Action rotary files allow you to move comfortably in the canal thanks to the special structure of its surface, while its variable taper is designed to evacuate debris from the canal. The developed special heat treatment technology ensures that the Ni-Ti wire has a much more flexible and durable structure and also turns it into a Gold or Blue color. EndoArt Action Gold files have a higher sharpness efficiency and higher fracture resistance than standard Ni-Ti files. EndoArt Action Blue files allows pre-bending before entering the canal thanks to its high flexibility efficiency with special control memory and it has a much higher fracture resistance than standard Ni-Ti files. It adapts to the shape of the canal in curvature structures.

Action Kit    SX - S1 - S2 - F1 - F2 - F3

**Action Gold Rotary Files**

- Rotation Files
- Variable Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

Control Memory  
Heat Treated Wire

Convex Triangle  
Cross - Section



**Action Blue Rotary Files**

- Rotation Files
- Variable Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

**Wire** : Gold  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files

**Wire** : Blue  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files



**WhatsApp**  
Get information!

**SMART**

04 Taper	Starter		Finishing					
	Torque: 3.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm

04 Taper	Starter		Finishing					
	Torque: 3.0 Ncm 300 Rpm	Torque: 1.6 Ncm 300 Rpm	Torque: 1.6 Ncm 300 Rpm	Torque: 1.6 Ncm 300 Rpm	Torque: 1.6 Ncm 300 Rpm	Torque: 1.6 Ncm 300 Rpm	Torque: 1.6 Ncm 300 Rpm	Torque: 1.6 Ncm 300 Rpm

EndoArt Smart rotary files have fixed %4 taper allow you to work more safely in narrow canals. The developed special heat treatment technology ensures that the Ni-Ti wire has a much more flexible and durable structure and also turns it into a Gold or Blue color. In addition to its high sharpness efficiency, EndoArt Smart Gold has a higher fracture resistance than standard Ni-Ti files. EndoArt Smart Blue allows pre-bending before entering the canal thanks to its high flexibility efficiency with special control memory and it has a much higher fracture resistance than standard Ni-Ti files. It adapts to the shape of the canal in cases with high slopes.

Smart Kit 15/04-20/04-25/04-30/04-35/04-40/04 19/08-10/04-15/04-20/04-25/04-30/04 19/08-15/04-20/04-25/04-30/04-35/04 19/08-16/02-20/04-25/04-30/04-35/04

**Smart Gold Rotary Files**

- Rotation Files
- Regular Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

Control Memory  
Heat Treated Wire

Convex Triangle  
Cross - Section



**Smart Blue Rotary Files**

- Rotation Files
- Regular Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

**Wire** : Gold  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files

**Wire** : Blue  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files



**WhatsApp**  
Get information!



TOUCH

10/06	15/06	20/06	25/06	30/06	35/06	40/06	50/06
Shaping		Finishing					
Torque: 3.0 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm	
Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm	

Excellent Cutting Efficiency  
Maximum File Flexibility

10/06	15/06	20/06	25/06	30/06	35/06	40/06	50/06
Shaping		Finishing					
Torque: 3.0 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm	
Torque: 1.8 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm	

EndoArt Touch rotary files fixed %6 taper allow you to prepare at a standard angle along the canal. The developed special heat treatment technology ensures that the Ni-Ti wire has a much more flexible and durable structure and also turns it into a Gold or Blue color. In addition to its high sharpness efficiency, EndoArt Touch Gold has a higher fracture resistance than standard Ni-Ti files. EndoArt Touch Blue allows pre-bending before entering the canal thanks to its high flexibility efficiency with special control memory and compared to standard Ni-Ti files, it has a much higher fracture resistance. It adapts to the shape of the canal in cases with high slopes.

Touch Kit 15/06-20/06-25/06-30/06-35/06-40/06 19/08-15/04-20/04-20/06-25/06-30/06 19/08-10/06-15/06-20/06-25/06-30/06 19/08-15/06-20/06-25/06-30/06-35/06

**Touch Gold Rotary Files**

- Rotation Files
- Regular Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box



**Touch Blue Rotary Files**

- Rotation Files
- Regular Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

**Wire** : Gold  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files

**Wire** : Blue  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files



WhatsApp  
Get information!

# ENDOART®



FINAL TOUCH



**Reciprocation**



EXPERT




EndoArt Expert files work with a reciprocating motion that allows you to perform root canal treatment with a single file. It includes alternative tapers, %6 and %4 angles for single-file systems, and allows you to prepare at a standard angle along the canal. EndoArt Expert files, which can be used with endomotors capable of reciprocating movement, work with a rotation angle adjustment of 150° left and 30° right. The developed special heat treatment technology makes the Ni-Ti wire a much more flexible and durable structure. EndoArt Expert Gold files have higher sharpness and higher fracture resistance than standard Ni-Ti files. EndoArt Expert Blue allows pre-bending before entering the canal thanks to its high flexibility efficiency with special control memory and compared to standard Ni-Ti files, it has a much higher fracture resistance.

Expert Kit R25-R25-R25-R25-R40-R50 R25-R25-R40-R40-R50-R50 R25-R25-R40-R50 R25-R40-R50

**Expert Gold Reciprocation Files**

- Reciprocation Files
- Regular Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box



**Wire** : Gold  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files

**Expert Blue Reciprocation Files**


- Reciprocation Files
- Regular Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

**Wire** : Blue  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files





PEDO



EndoArt Pedo rotary files are produced specially for pediatric cases and has a length of 18 mm. The improved heat treatment technology of EndoArt Pedo ensures that the Ni-Ti wire has a much more flexible and durable structure. EndoArt Pedo files have higher sharpness and fracture resistance than standard Ni-Ti files.

Pedo Kit 15/06-20/04-25/04 15/06-20/04-30/04

**Pedo Gold Rotary Files**

- Rotation Files
- Specially produced for pedodontic cases
- Regular Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

Control Memory  
Heat Treated Wire

Convex Triangle  
Cross - Section



**Pedo Blue Rotary Files**

- Rotation Files
- Specially produced for pedodontic cases
- Regular Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

**Wire** : Gold  
**Length** : 18mm  
**Packaging** : 3 Files

**Wire** : Blue  
**Length** : 18mm  
**Packaging** : 3 Files



WhatsApp  
Get information!

## RETREATMENT



EndoArt Retreatment rotary files developed for repeated root canal treatments. Thanks to its high cutting efficiency, it is used in the process of removing the materials in the canals. The system consists of three instruments, D1 - D2 and D3 which are designed for cleaning the canal filling in the (D1) coronal, (D2) middle, and (D3) apical thirds. While D1 has an active tip at the first entry into the root canal, the tips of continuing instruments are not cutting.

Retreatment Kit D1-D1-D2-D2-D3-D3

### Retreatment Files

Rotation Files  
 Specially produced for retreatment cases  
 Designed to remove filling materials like Gutta Percha  
 x2 times better Cyclic Fatigue Resistance than conventional Files  
 Better Cutting Efficiency  
 Great Fracture Resistance  
 Lasered Taper and Size information on the shank  
 Autoclavable Box

**Wire** : Gold  
**Length** : 16mm, 18mm, 22mm  
**Packaging** : 6 Files



## PATH FILE




Path files are used to create a glide (Glide) path. The use of angled files in combination with hand files before use, reduces the risk of step formation and perforation. The size and taper of the file are 16/02. The ends of the file are rounded to prevent the formation of step and zip formation.

### Path Gold Rotary Files

Rotation Files  
 Specially produced for creating glide path  
 Easily reach apex  
 Better Flexibility  
 x2 times better Cyclic Fatigue Resistance than conventional Files  
 Great Fracture Resistance  
 Lasered Taper and Size information on the shank  
 Autoclavable Box

**Wire** : Gold  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files



WhatsApp  
 Get information!





## EndoArt K Files

Codes : ■06 ■08 ■10 ■15 ■20 ■25 ■30 ■35 ■40 □45  
■50 ■55 ■60 ■70 ■80 □90 ■100 ■110 ■120 ■130 ■140

Ass. 15/40 – Ass. 45/80 – Ass. 90/140

Stainless Steel

Length: 21 mm, 25 mm, 31 mm

6 pieces in a box



## EndoArt H Files

Codes : ■06 ■08 ■10 ■15 ■20 ■25 ■30 ■35 ■40 □45  
■50 ■55 ■60 ■70 ■80 □90 ■100 ■110 ■120 ■130 ■140

Ass. 15/40 – Ass. 45/80 – Ass. 90/140

Stainless Steel

Length: 21 mm, 25 mm, 31 mm

6 pieces in a box



## EndoArt Spreader

Codes: ■15 ■20 ■25 ■30 ■35 ■40/ Ass. 15/40

Stainless Steel

Length: 21 mm, 25 mm

6 pieces in a box



## EndoArt C Files

Codes: ■06 ■08 ■10 ■12.5 □15/ Ass. 06-10

Stainless Steel

Length: 25 mm

6 pieces in a blister package

## EndoArt Gutta Percha Points

Gutta Percha root canal filling material for easy fillings. It has variable taper and sizes.

**Properties**

- 02/04/06 Taper Options
- Excellent Radiopacity
- Length Marked
- Hand Rolled



## EndoArt Paper Points

Absorbent Paper Point to dry root canals.

**Properties**

- Highly Absorbent
- Sterile Packaging
- Variety Sizes and Tapers
- Length Marked
- Hand Rolled





## EndoArt Fiber Post

Universal resin glass fiber post system for esthetic and metal free restorations. It has different type of diameters which are 1,2 mm, 1,4 mm, 1,6 mm and 1,8 mm.

### Properties

Radiopaque  
 4 drills available inside of the kit  
 Conical tip  
 High resistance  
 Colour coded  
 It is compatible with adhesive resin cements

### Packaging

EndoArt Fiber Post Kit – 5 pcs. 1,2 mm posts, 5 pcs. 1,4 mm posts, 5 pcs. 1,6 mm posts, 5 pcs. 1,6 mm posts, 4 drills assorted  
 Ø 1.2 - Ø 1.4 - Ø 1.6 - Ø 1.8

EndoArt Fiber Post Refill - 5 pcs. 1,2 mm posts

EndoArt Fiber Post Refill - 5 pcs. 1,4 mm posts

EndoArt Fiber Post Refill - 5 pcs. 1,6 mm posts

EndoArt Fiber Post Refill - 5 pcs. 1,8 mm posts



## RubyCal Calcium Hydroxide Paste With Barium Sulphate

Calcium hydroxide premixed paste non-setting for root canal therapy and vital pulp therapy. It comes with application tips.

### Properties

- High Radiopacity
- Easy removal
- High pH > 12
- Excellent handling characteristics

### Packaging

- 2 gr syringe
- 2 gr syringe x 2
- 2 gr syringe x 4



## RubyCal Calcium Hydroxide Paste With Iodoform

Calcium hydroxide premixed paste non-setting for root canal therapy and vital pulp therapy. It comes with application tips.

### Properties

- Radiopacity
- Promote apexification
- Disinfects Canal
- Excellent handling characteristics

### Packaging

- 2 gr syringe
- 2 gr syringe x 2
- 2 gr syringe x 4



**ruby dent**





incidental



## RubyFill Nano Composite

RubyFill Nano is radiopaque nano hybrid composite system which can be the photocurable applied in the anterior and posterior field.

### Properties

- Universal Nano Composite
- It has light refraction property like a normal tooth due to the nano fillers in it
- It has high filler rate
- It can be polished high with high degree
- It has high resistance to abrasion
- It provides ease of use

### Packaging

- Refill – 4 gr in a syringe
- Kit 1 - 4 gr syringe x 6 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories
- Shades – A1 – A2 – A3 – A3,5 – B1 – B2 – C1 – D3 – OA2 – Extra Bleach



## RubyComp Nano Composite

RubyComp Nano is radiopaque nano hybrid composite system which can be the photocurable applied in the anterior and posterior field.

### Properties

- Higher Esthetic
- Low shrinkage
- It has high filler rate
- It can be polished with high degree
- It has high resistance to abrasion
- It provides ease of use
- Long durability
- Non sticky

### Packaging

- Refill – 4 gr in a syringe
- Kit 1 - 4 gr syringe x 6 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories
- Shades – A1 – A2 – A3 – A3,5 – B1 – B2 – C1 – D3 – OA2 – Extra Bleach



ISO 4049

CLASS 2A



WhatsApp  
Get information!



## RubyComp P Posterior Composite

RubyComp P is micro hybrid composite system which can be the photocurable applied in the posterior field.

### Properties

- Non stick to instruments
- It has surface hardness is like real enamel
- It has high intensity
- It has a perfect consistence can be shaped/formed easily
- It has high filler rate
- It resists to chewing force
- There is not discoloration after the polimeration.

### Packaging

Refill – 4 gr in a syringe  
 Kit 1 - 4 gr syringe x 6 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories  
 Shades – A1 – A2 – A3 – A3,5



## RubyFlow Flowable Composite

RubyFlow is flowable light curing radiopaque composite.

### Properties

- Low shrinkage, low wear
- High gloss polishability
- Beautiful esthetic restorations
- Very low surface roughness
- Excellent handling
- Non sticky
- Very good marginal adaptation

### Packaging

Refill – 2 gr in a syringe  
 Kit 1 - 2 gr in a syringe x 2 pcs.  
 Kit 2 - 2 gr syringe x 4 pcs.  
 Kit 3 - 2 gr syringe x 4 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories  
 Shades – A1 – A2 – A3 – A3,5 – B1 - OA2



ISO 4049

CLASS 2A



WhatsApp  
Get information!





## RubyBond Adhesive System

RubyBond is light curing single component 5. generation bonding agent for the adhesive restorations. It universal and compatible with all composites.

### Properties

Low abrasion  
Elastic connection  
5. generation

### Packaging

5 ml in a bottle  
5 ml in a bottle + 3ml Etching



## RubySEBond Self Etching Adhesive System

RubySEBond is self-etching light curing one component 7 generation adhesive for the adhesive restorations.

### Properties

Very strong bonding to enamel and dentine  
No extra etching with phosphoric acid  
Long lasting restorations  
Time savings  
7. generation

### Packaging

5 ml in a bottle



CLASS 2A



 **WhatsApp**  
Get information!



## RubyEtch Etching Gel

RubyEtch is an etching gel containing 37% phosphoric acid. Delivered in 3 ml and 50 ml syringe with very thin application needles.

### Properties

Easily Rinse  
Blue coloured  
Perfect Surface Retention

### Packaging

3 ml syringe  
3 ml syringe x 2  
3 ml syringe x 4  
50 ml syringe



CLASS 2A



 **WhatsApp**  
Get information!



## RubyLiner with Glass Ionomer

Light curing glass ionomer based cavity liner to protect pulp of the tooth. It is used under direct and indirect restorations, like composites, amalgams, ceramics and another metal restorations.

### Properties

Excellent biocompatibility  
Releases flouride  
Radiopaque  
Excellent handling characteristics

### Packaging

2 gr syringe  
2 gr syringe x 2  
2 gr syringe x 4



## RubyCal LC Calcium Hydroxide

Light curing calcium hydroxide for pulp capping and liner material under dental filling.

### Properties

Direct & indirect pulp capping  
Light curing  
Single component

### Packaging

2 gr syringe  
2 gr syringe x 2  
2 gr syringe x 4



ISO 4049

CLASS 2A



WhatsApp  
Get information!



## Ruby Pit and Fissure Sealant

Light curing sealant with fluoride releasing for occlusal surfaces. It can be used as preventive products.

### Properties

Fluoride releasing  
Radiopaque  
Easy application  
Low shrinkage

### Packaging

2 gr syringe  
2 gr syringe x 2  
2 gr syringe x 4



## FlexDam Gingiva Barrier

Light curing resin-based gingiva barrier to protect soft tissue.

### Properties

Easily removable after treatment  
It minimizes heat and soft tissue sensitivity during curing  
Blue shade  
5 second curing

### Packaging

2 gr syringe  
2 gr syringe x 2  
2 gr syringe x 4



CLASS 2A



WhatsApp

Get information!




## RubyShine Prophylaxis Paste

RubyShine Prophylaxis Paste is specially formulated of polishing and cleaning agents which has nano particuls in it. Available in mint, strawberry, bubblegum and green apple flavor.

### Packaging

30 gr in a bottle  
150 gr in a bottle  
340 gr in a bottle



Green Apple

Mint

Bubblegum

Strawberry



Class I



WhatsApp  
Get information!

## RubyPlaton Polishing Discs

Has different kind of discs which has different grins and different sizes. It allows excellent results in polishing. Flexible and easy to use


### Properties

14 mm blue, yellow, white, green, brown, pink

- **White** – Extra Fine for Polishing
- **Yellow** – Fine for Finishing
- **Green** – Medium for Countoring
- **Blue** – Coarse for Gross Reduction
- **Brown** – X Coarse for Gross Reduction
- **Pink** – XX Coarse for Gross Reduction

### Packaging

80 pcs discs and 1 mandren in a plastic box  
240 pcs discs and 4 mandren in a plastic box



## RubyPlaton Polishing Strips

Has different kind of strips which has different grin sizes. They are produced combining two different types of grins.

### Properties

It allows excellent results in polishing. Flexible and easy to use

- **White** – Extra Fine for Polishing
- **Yellow** – Fine for Finishing
- **Green** – Medium for Countoring
- **Blue** – Coarse for Gross Reduction

### Packaging

75 pcs strips in a plastic box  
80 pcs discs and 1 mandren + 75 pcs strips in a plastic box  
240 pcs discs and 4 mandren + 75 pcs strips in a plastic box





## RubyPlaton Wooden Wedges

RubyPlaton wooden wedge it has anatomic form for you to make excellent filling.

### Properties

- Superthin short
- Thin short
- Thin long
- Medium short
- Medium long

### Packaging

100 pcs. wedges in a plastic box  
 200 pcs. wedges in a plastic box  
 500 pcs. wedges in a plastic box



## RubyPlaton Transparent Wedges

RubyPlaton transparent wedge is formed by special transparent polymer for excellent filling.

### Properties

Not easy to break long lived with high quality and durability  
 Easy to use  
 Anatomic form  
 Small size - Medium size

### Packaging

100 pcs in a plastic box

CE Class I



WhatsApp  
 Get information!



## RubyMark Articulating Paper Straight

Great sensitive articulating paper to get simple marking. It is straight form on the 40µ, 80µ and 200µ.

### Properties

Non sticky  
Smooth surface  
Double sided options.

### Packaging

40µ – Blue/Blue Straight 144 pcs. paper  
40µ – Red/Red Straight 144 pcs. paper  
40µ – Blue/Red Straight 144 pcs. paper  
80µ – Blue/Blue Straight 144 pcs. paper  
80µ – Red/Red Straight 144 pcs. paper  
80µ – Blue/Red Straight 72 pcs. paper  
80µ – Blue/Red Straight 144 pcs. paper



Class I



WhatsApp

Get information!





## RubyMark Articulating Paper Horseshoe

Great sensitive articulating paper to get simple marking. It is horseshoe form on the 80µ.

### Properties

Non sticky  
Smooth surface  
Double sided options.

### Packaging

80µ – Blue/Red Horseshoe 72 pcs. paper

## RubyMark Articulating Paper 200 Micron Straight

Great sensitive articulating paper with special fleece paper to get simple marking. It is straight form on the 200µ.

### Properties

Non sticky  
Smooth surface  
Most Accurate Occlusion

### Packaging

200µ – Blue Straight 300 pcs. paper  
200µ – Red Straight 300 pcs. paper



CE Class I



WhatsApp  
Get information!



## RubyCrownRemover

It is practical crown remover for easy procedure. It removes crown using the way of impact method with arc.

### Properties

- Can be sterilized at 134 °C
- 3 types of remover tips
- Different types of impact



## RubyPlaton Cotton Roll

High quality cotton rolls made from %100 cotton. They are non-linting, non-sterile, soft texture rolls.

### Properties

- It retracts cheek effective
- It does not adhere to mucous membrane
- Highly absorbent
- Size: 4 x 1 cm
- Maintain shape when wet

### Packaging

- No: 1 / 500 pieces bag
- No: 2 / 500 pieces bag
- No: 3 / 500 pieces bag

CE Class I



WhatsApp  
Get information!



## Arma Dental Study Model


Dental study model for students. It shows great and real-like experience during application.

### Properties

Hard Plastic jaw base  
 Non-breakable  
 The students thus practice under "human" conditions  
 Custom made screws for tooth fixing  
 Natural anatomy and occlusion

### Packaging

1 piece dental study model with screwed teeth



## Arma Dental Phantom Teeth

Dental study model teeth for students. It shows great and real-like experience during application. Teeth are made from melamine material.

### Properties

Custom made screws for tooth fixing  
 Natural anatomy and occlusion

### Packaging

1 set dental study model teeth and screws with teeth chart

# ENDOART®

FINAL TOUCH



incidental